

» Passion
for **Precision** «



MPS

Multi Platform Systems

IMPLANT

PICK-UP

Closed tray



»» Welcome to MEDENTIKA® ««

Anyone who meets MEDENTIKA® personnel will get to know a team. A team that since the foundation of the company collectively contributes knowledge and experience. A team that perfects products with élan. On a daily basis and with passion.

MEDENTIKA® has character. And it is this character that gives us a morally sound base. In a world that is becoming more and more focused on the laws that drive market forces, we feel that this position is not only right, but important. What characterises MEDENTIKA® is also our belief in the power of our ideas and in the innovative power of our products. We are MEDENTIKA® – and if the MEDENTIKA® label is on the packet, there is more than just an implant component inside, but an implant component that bears the knowledge and the experience that we have been putting into our products since MEDENTIKA® was founded. We offer our customers our enthusiasm – a feature we use to perfect our products on a daily basis.

We take this task seriously, but we also embrace it with passion. And the result is an original product manufactured with the highest interface precision from an extremely wide range of prosthetics.

Our catalogues reflect our attitude in every single product.

Our prosthetic components offer you new ways and even more treatment opportunities providing your patients more comfort and functionality!

Sincerely your
Jürgen Kälber » Bernd Gaddum » Frank Fix
and the entire MEDENTIKA® team.

**»» implant pick-ups
for closed
impression ««**

NEU

COLOUR CODE

The implant pick-up for closed impression are colour-coded corresponding to the colour of the respective implant manufacturer.

POSITIONING CAP

The implant pick-ups are delivered including the positioning cap, which guarantees the exact and precise repositioning of the implant pick-up into the impression.

SERIES ASSIGNMENT

For better visibility the implant pick-ups are additionally marked with the serial letter of the respective MEDENTIKA® series.




IMPLANT CONNECTION

The implant pick-up is made of titanium Grade 5 CF and has a highly precise implant connection. The implant pick-up will be delivered including positioning cap and retention screw.

>> E-Series <<< Nobel Biocare NobelReplace Tapered

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF



Implant connection	NP 3,5	RP 4,3	WP 5,0
Article-No.	E 30	E 31	E 32

>> H-Series <<< BIOMET 3i Certain

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF



Implant connection	D 3,4	D 4,1	D 5,0
Article-No.	H 30	H 31	H 32

>> F-Series <<< Nobel Biocare NobelActive - NobelReplace Conical

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF



Implant connection	D 3,0	NP 3,5	RP 4,3/5,0
Article-No.	F 32	F 30	F 31

>> L-Series <<< Straumann Bone Level

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF



Implant connection	NC 3,3	RC 4,1/4,8
Article-No.	L 30	L 31

>> N-Series << Straumann Tissue Level

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF




Implant connection	NNC 3,5	RN 4,8	WN 6,5
Article-No.	N 33	N 31	N 32

>> S-Series << DENTSPLY Implants ASTRA TECH OsseoSpeed TX

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF




Implant connection	D 3,5/4,0	D 4,5/5,0
Article-No.	S 30	S 32

>> R-Series << Zimmer Dental Tapered Screw-Vent MIS SEVEN Internal Hex BioHorizons Tapered Internal/Internal Plus/Tissue Level

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF




Implant connection	D 3,5	D 4,5	D 5,7
Article-No.	R 30	R 31	R 32

>> T-Series << DENTSPLY Implants XIVE S

Implant pick-up
closed tray

- incl. retention screw
- incl. positioning cap
- Titanium Grade 5 CF



Implant connection	D 3,4	D 3,8	D 4,5	D 5,5
Article-No.	T 30	T 35	T 31	T 32

IMPRINT

Editor: MEDENTIKA® GmbH
Hammweg 8-10
76549 Hügelshiem
+49 (0)7229 69912-0
+49 (0)7229 69912-20
info@medentika.de
www.medentika.de

Design: Der WeberFink GbR
Graphic Design Studio
www.weberfink.de

MEDENTIKA® GmbH
Hammweg 8-10
76549 Hügelshheim
Tel.: +49 (0)72 29 69912-0
Fax: +49 (0)72 29 69912-20
info@medentika.de
www.medentika.de

»» Passion
for **Precision** ««

